---

# IGN<sup>®</sup> Calbénium<sup>®</sup>

---

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*IGN EVO® is a water decontamination system fully automatic and adaptable to all dental units. The water to be treated is enriched with CALBENIUM® thanks to an automatic system.*

*It is used with CALBENIUM®, a decontamination product made of ultrapure components.*


*The CALBENIUM® is used in liquid form. Decontamination is obtained through ionization, created by a magnetic field and a silver induction.*

*This system insure also a disconnection of the water supply with the water sprays, in conformity with norm NF EN 1717 (protection against the drinking water pollution). This model is provided with a regulator, two filters and external quick connections. Alarm systems are controlling the CALBENIUM® level.*

*IGN EVO® is connected to water, air and power supplies.*

*Design, IGN EVO® is perfectly adaptable in a dental office. IGN EVO® is inside a housing and may be fitted easily next to or onto the dental units.*

*Every 2 weeks, an using with a CALBENIUM® dose multiplied by 5 produce a "bacteriological reset". It allows a correct maintenance conform to ISO norms.*

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*Composition : EDTA - benzalkonium chloride - sodium tosyl-chloramide - allantoin - aspartam - spearmint or lemon - sorbitol . CALBENIUM® is made with "ultrapur" products containing less than 0,1 % of impurities.*

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*SPRAY MIST : A DROPLET CONTAMINATION RISK*  
*Various risk of contamination are present in a dentist's surgery : from the dental care unit, from patient to patient, from patient to dentist, from instruments to the assistant.*

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*Created in 1988, it was the first of the range, a semi-automatic one. The ionization process requires a chrome electrode.*

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*Created in 1993, it is a fully automatic device. The ionization is produced by an electronic PC Board and a silver electrode.*

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*Created in 1999, the ionization in this version is obtained by a silver electrode and a permanent battery of 1.5 V.*

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*Created in 2011, it is a fully automatic device. The ionization is obtained through a magnetic field and a silver induction.*

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